

# HASO SWIR R.FLEX2

**FOR  $\lambda = 1.0 - 1.7 \mu\text{m}$   
WITH  $\lambda/200$  RMS ACCURACY**

**LARGE RANGE OF F/#  
AVAILABLE**

**COMPACT AND ROBUST  
FOR EASY INTEGRATION**

**UP TO 150 Hz  
ACQUISITION FREQUENCY**




SWIR wavefront and MTF measurements for characterizing optical components, such as lenses, filters, waveplates, telescopes and complex optical systems

## A UNIQUE SET OF ADVANTAGES

- $\lambda/200$  rms measurement accuracy in double-pass configuration
- Collimated or diverging exit beam with the F/# that matches to optics under test
- Highly accurate wavefront analysis even with central obscuration and/or spider-beam types
- Insensitive to vibrations and atmospheric turbulences
- Removable wavefront sensor for using it as a stand-alone unit
- Several accessories available, such as laser diode light sources, reference mirrors for calibration and translation stages


Example of measurement configuration



光貿易株式会社 〒113-0034 東京都文京区湯島3-13-8湯島不二ビル301  
メール: [contact@hikari-trading.com](mailto:contact@hikari-trading.com) ウェブサイト: <http://www.hikari-trading.com>  
TEL: 03-3832-3117 FAX: 03-3832-3118


### Measuring large concave mirrors

HASO R-Flex2 SWIR has been optimized using proprietary designs that enable manufacturers to accurately measure large uncoated concave mirrors by positioning the unit to measure at the center of curvature.




### Characterizing complex optical systems

Complex optical systems such as telescopes and collimators can be readily characterized by HASO R-Flex2 SWIR. The best focal point can be found using wavefront error whereas, if the focus point is defined mechanically, optics can be aligned for that point.




### Measuring lenses on-axis

Any diameter lenses are easily measured with HASO R-Flex2 SWIR by using a coated or uncoated flat reference mirror to reflect the beam back to the wavefront sensor without adding any aberrations.




### Characterizing lenses in the field

By mounting the HASO R-Flex2 SWIR onto a translation stage and orienting the flat reference mirror correspondingly, you can qualify lenses at any point in the field.




### Characterizing & aligning beam expanders

HASO R-Flex2 SWIR 's modularity is particularly useful since its focusing unit dismounts quickly and therefore a collimated beam can be used as an illumination source to characterize the beam expander without adding any aberrations.



### Working with external sources

High N/A external sources can be accurately measured because the optical head can be completely characterized (left image). Dismount it, and you can use the wavefront sensor as a stand-alone unit (right image).



HASO SWIR



Wavefront sensor

R-FLEX2



Collimator

MOD



Focusing module



Optics under test



Reference flat mirror

## R-Flex2 SWIR SPECIFICATIONS


Compatible wavefront sensor	HASO SWIR
Wavefront measurement accuracy in double pass configuration	$\lambda/200$ RMS
Aperture dimension	9.30 x 7.44 mm <sup>2</sup>
Number of phase points	40 x 32
Collimated beam diameter	~12mm*
Maximum acquisition frequency (Hz)	150 Hz

## FOCUSING MODULES

Focusing module name	F number	Focal length (mm)	Wavefront error (nm RMS)**	Working distance from the module (mm)	Required back power (%)
MOD F20	2.7	20	250	10.0	3
MOD F31	4.2	31	150	10.0	3
MOD F40	5.4	40	100	3.3	3
MOD F50	6.8	50	25	3.3	3
MOD F60	8.1	60	50	12.8	3
MOD F76	10.2	75.3	25	12.8	3

\* Output light beam diameter only, its useable area is always determined by sensor's aperture dimension.

\*\* WaveFront Error (WFE) at the output of the module for a circular pupil corresponding to the nominal F/#



Screenshot of WaveView software showing reconstructed wavefront of a beam with an obscuration.

Features	Values	Name	Values (um)
Maximize pupil	Zernike	Z01 Tilt 0°	0.0000
Coefficients type	22	Z02 Tilt 90°	0.0000
Number of coefficients	RMS	Z03 Focus	-0.1000
Coefficients normalization	0.0000	Z04 Astig. 0°	0.0200
Pupil X center (mm)	7.304	Z05 Astig. 45°	-0.0100
Pupil Y center (mm)		Z06 Coma 0°	0.0000
Pupil radius (mm)		Z07 Coma 90°	0.0000
		Z08 Spherical	0.0000
		Z09 Trefoil 0°	0.0000
		Z10 Trefoil 90°	0.0000
		Z11 5th Astig. 0°	0.0400
		Z12 5th Astig. 45°	0.0000
		Z13 5th Coma 0°	0.0000
		Z14 5th Coma 90°	0.0000
		Z15 5th Spherical	0.0000
		Z16 Trefoil 0°	0.0000
		Z17 Trefoil 90°	0.0000

Zernike coefficient window

## Accessories

### Translation stages

Our  $\Theta X\Theta Y$  rotation stage for angular alignment or the 5-axis stage that provides 2-way rotation around X and Y axes as well as 3-way translation along X, Y and Z axes is a perfect complement to the HASO R-Flex system.

### Software add-on

HASO R-Flex is delivered with WaveView software, which is a leading wavefront metrology software providing 180 independent features. We also offer optional software modules including MTF (Modulation Transfer Function) and PSF (Point Spread Function) that increase the functionality of HASO R-Flex system.

### Reference mirrors

Spherical reference mirror ( $\varnothing 20\text{mm}$  useful pupil,  $R=15\text{mm}$ ,  $F/0.75$ ) for the calibration of HASO R-Flex in double-pass measurement configuration

Flat reference mirror for autocollimation. Several options are available in diameter and flatness.

### Single-Mode Laser Source (SMLS)

For those who want to use their HASO R-Flex at different wavelengths, we provide additional single-mode diode lasers to further expand the versatility of the system. Please contact us if you prefer to use your own light source.

### NEW: R-Flex Kit

Kit for R-Flex calibration check and maintenance. It includes a light source for pre-alignment and fiber checking, a dust remover, a telescopic mirror, a torch lamp, a plane mirror, a retroreflector, and the instruction for R-Flex calibration check.

### Available SMLS wavelengths:

Model name	Wavelength (nm)	Maximal power (mW)
SMLS 1064-S	1064	4.5
SMLS 1550-S	1550	4.5
SMLS custom	Ask	Ask



SMLS

### SWIR related products

- HASO SWIR
- HASO4 SWIR 1550
- R-Flex2 SWIR
- R-Flex Large Aperture (LA) SWIR



光貿易株式会社 〒113-0034 東京都文京区湯島3-13-8湯島不二ビル301  
メール: [contact@hikari-trading.com](mailto:contact@hikari-trading.com) ウェブサイト: <http://www.hikari-trading.com>  
TEL: 03-3832-3117 FAX: 03-3832-3118